

Akriti




Surgical Products

About Us: Akriti is a group of companies Based in Hyderabad and engage in Eye Care products development, manufacturing and trading facility across India, Asia, Africa, United Kingdom and United States, having its own office in various Locations Headed by Dr. Kuldeep Raizada, Ph D, BCO, BADO, FAAO, Who is pioneer in Prosthetics Graduated from USA, Akriti is located in medical devices park, In Hyderabad spread over 75000 Sft.

Vision: Manufacturing facility to develop & Manufacturer the quality products for the develop & developing world, which can be affordable for everyone.

Statement: Ensure to timely deliver a quality products to everyone.

Our Manufacturing Facility at Hyderabad



Our Quality Certifications

Certificate

This is to Certify that WHO-Good Manufacturing Practice of the

AKRITI OPHTHALMIC PRIVATE LIMITED

M34, M35, M36, Medical Device Park, Sultanpur Village, Pathancheru Mandal, Hyderabad, Telangana, 502319, India

Has been assessed and found to be in accordance with the requirements standard detailed below

WHO-GMP

This Certificate is valid for the following scope

Manufacture of Medical Devices, Surgical Instruments, Ophthalmic Consumables & Hospital Furniture.

Certificate No : WHO-MMNXXII-04-11484
 NACE : DL-33-19
 Issue Date : 05-APR-2022
 Valid Until : 04-APR-2025

Adpanghal
 Authorized Signatory
 WRG Certification

This certificate remains the property of American International Accreditation-Bureau of Accredited Register, (USA)

Certificate

This is to Certify that Medical Device Quality Mangement System of the

AKRITI OPHTHALMIC PRIVATE LIMITED

M34, M35, M36, Medical Device Park, Sultanpur Village, Pathancheru Mandal, Hyderabad, Telangana, 502319, India

Has been assessed and found to be in accordance with the requirements standard detailed below

13485:2016

This Certificate is valid for the following scope

"Manufacture of Ophthalmic Medical Devices, Surgical Instruments, Ophthalmic Consumables, Disposable, Eye Wear, Safety wear & Hospital Furniture"

Certificate No : MMN-MMNXXII-04-11485
 NACE : DL-33-19
 Issue Date : 05-APR-2022
 Valid Until : 04-APR-2025

Adpanghal
 Authorized Signatory
 WRG Certification

This certificate remains the property of American International Accreditation-Bureau of Accredited Register, (USA)

51021006051

WRG CERTIFICATIONS

Certificate of Registration

WRG Certifications hereby certifies that the Quality Management System of:

AKRITI OPHTHALMIC PRIVATE LIMITED

M34, M35, M36, Medical Device Park, Sultanpur Village, Pathancheru Mandal, Hyderabad, Telangana - 502319, India

Has been assessed and found to operate in compliance and meets the requirement of following standard

ISO 9001:2015

for the scope of

Manufacture, Trade and Supplier of Medical Devices, Surgical Instruments, Ophthalmic Consumables & Hospital Furniture.

Valid date of certification : Apr 04, 2022
 Current date of certification : Apr 04, 2022
 Date of expiry : Apr 04, 2025

1st Surveillance Date : Mar 2023
 2nd Surveillance Date : Mar 2024

Certificate Number : GMS/M00004-11485
 NACE Code : DL33-10

WRG Certification
 Adpanghal
 Authorized Signatory

AIAO-BAR
 QMS
 ACCREDITED
 USA

WRG
 ISO 9001
 CERTIFIED

This certificate remains the property of American International Accreditation-Bureau of Accredited Register, (USA) & WRG
 Register Group (certification) shall be withdrawn. This certificate is valid only from the date of the issue. Near to renewal, it will be re-issued if compliance is maintained. Compliance & Accreditation by AIAO-BAR. www.aiao-bar.com Registration No - AIAO-BAR-0000033-0



Certificate

This is to Certify that WHO-Good Manufacturing Practice of the

AKRITI OPHTHALMIC PRIVATE LIMITED

M34, M35, M36, Medical Device Park, Sultanpur Village, Pathancheru Mandal,
Hyderabad, Telangana, 502319, India

*Has been assessed and found to be in accordance with the
requirements standard detailed below*

WHO-GMP

This Certificate is valid for the following scope

Manufacture of Medical Devices, Surgical Instruments,
Ophthalmic Consumables & Hospital Furniture.

Certificate No : WHO-MMXXII-04-11484
NACE : DL-33.10
Issue Date : 05-APR-2022
Valid Until : 04-APR-2025

A.dpanghal
Authorised Signatory
WRG Certification

This certificate remains the property of American International Accreditation Organization-Bureau of Accredited Registrars, (USA)



Validity of the certificate is subject to successful completion of surveillance audits. Further clarification regarding the scope of this certificate and the applicability of standard may be obtained by consulting the WRG Certification Organization, India Aditya Corporate Hub, Raj Nagar, Ghaziabad-201002 (U.P.) India, E-mail : wrgcert@gmail.com, website : www.wrgcert.com, www.aiao-bar.org



Certificate

This is to Certify that Medical Device Quality Mangement System of the

AKRITI OPHTHALMIC PRIVATE LIMITED

M34, M35, M36, Medical Device Park, Sultanpur Village, Pathancheru
Mandal, Hyderabad, Telangana, 502319, India

*Has been assessed and found to be in accordance with the
requirements standard detailed below*

13485:2016

This Certificate is valid for the following scope

"Manufacture of Ophthalmic Medical Devices, Surgical Instruments, Ophthalmic Consumables, Disposable, Eye Wear, Safety wear & Hospital Furniture"

Certificate No :MMS-MMXXII-04-11485
NACE :DL-33.10
Issue Date :05-APR-2022
Valid Until :04-APR-2025

A.dpanghal
Authorised Signatory
WRG Certification

This certificate remains the property of American International Accreditation Organization-Bureau of Accredited Registrars, (USA)



Validity of the certificate is subject to successful completion of surveillance audits. Further clarification regarding the scope of this certificate and the applicability of standard may be obtained by consulting the WRG Certification Organization, India Aditya Corporate Hub, Raj Nagar, ghaziabad-201002 (U.P.) India, E-mail : wrgcert@gmail.com, website : www.wrgcert.com, www.aiao-bar.org



Certificate Of Registration

WRG Certifications hereby certifies that the Quality Management System of:

AKRITI OPHTHALMIC PRIVATE LIMITED

*M34, M35, M36, Medical Device Park, Sultanpur Village, Pathancheru Mandal, Hyderabad,
Telangana - 502319, India*

has been assessed and found to operate in compliance and meets the requirement of following standard

ISO 9001:2015

for the scope of

Manufacture, Trader and Supplier of Medical Devices, Surgical Instruments, Ophthalmic Consumables & Hospital Furniture.

Initial date of certification : Apr 05, 2022
Current date of certification : Apr 05, 2022
Date of expiry : Apr 04, 2025

(Subject to surveillance mark present)

1st Surveillance Date : Mar 2023



2nd Surveillance Date : Mar 2024



Certificate Number : QMS-MMXXII-04-11483
NACE Code : DL-33.10

WRG Certifications

A.dpanghal

Authorised Signatory



Registry Information can be found at: www.wrgcert.com/register



ISO 9001

This certificate remains the property of American International Accreditation Organization – Bureau of Accredited Registrars, (AIAO-BAR) & World Registrar Group Certification, (WRG Certifications). This certificate is valid for three years from the date of its issuance. Must be returned, if certificate is withdrawn. Compliance & Accreditation by AIAO-BAR. www.aiao-bar.org. Registration No – AIAO-BAR-0305131-1

Forceps

K5-1560



Polack Double Corneal
Fixation Forceps

S0040



Castroviejo Forceps
0.12mm

S0141



Nugent Utility Forceps,
Serrated

S0041



Castroviejo Forceps
0.30 mm

S0021



Berke Ptosis Forceps

S0052



Colibri Forceps
1x2 teeth
0.40mm

S0030



Botvin Iris Forceps

S0053



Colibri Forceps
0.12 mm
Teeth

CE

6

Forceps

S0088



Harms-Colibri Forceps
0.12mm

S0078



Francis Chalazion Forceps

S0061



Dodick IOL Folding Forceps

S0095



Hunt Chalazion Forceps

S0065



Ernest Nucleus Cracker

S0084



S0084 Halsted Mosquito Forceps, Curved
S0085 Halsted Mosquito Forceps, Straight

S0066



Ernest-McDonald IOL Insertor

S0054



Conjunctiva Forceps

CE

7

Forceps

S0113



Lester Fixation Forceps

S0170 X



Tissue Forceps

S0171



Tissue Forceps
Delicate
Straight

S0132



McPherson Tying Forceps

S0168



Tennant Tying Forceps
Curved

S0177



Utrata Capsulorhexis Forceps

S0169



Tennant Tying Forceps
Straight

S0204



Sauer Fixation Forceps



Forceps

S0218



Moody Fixation Forceps
S0217 Left
S0218, Right

S0159



Squeeze Handle Forceps
S0159 serrated, straight
S0160 scissors, horizontal,
straight

S0226



Lims Forceps

S0204



Toothed Forcep

S0237



Hoskin forceps

S0170



Non Toothed Forceps

S0233



Superior Rectus Forceps

S0279



Cheatle Forceps



Forceps / Scissors

S0026



Bishop-Harmon Forceps

S0063



Dressing Forcep

S0162



Stern-Gills Capsulotomy Scissors

S0134



McPherson-Westcott Stitch Scissors

S0163



Stern-Gills Scissors

S0179



Vannus Scissors

Troutman-Katzin Scissors

S0174




S0174 Troutman-Katzin Scissors, small, left

S0175 Troutman-Katzin Scissors, small, Right



Scissors

S0016



Barraquer Iris Scissors

S0186



Westcott Tenotomy Scissors

S0164



Stevens Scissors
Curved

S0174



Troutman-Katzin Scissors
S0174 small, left
S0175 small, right

S0165




Stevens Scissors
Straight

585027



End Grip Scissor
585027 20 G
585028 23G

S0067



S0064 Enucleation Scissors, Cvd
0067 Eye Scissors, Cvd
S0068 Eye Scissors, Small
S0069 Eye Scissors, Straight
S0070 Eye Scissors, straight



11

Speculum

S0019



Barraquer Wire Speculum
S0019 Large
S0020 Small

S0115



Lester-Burch Eye Speculum

S0224




Muller Speculum

S0119



Lieberman Speculum
Solid blade

S0117



Lieberman Aspirating Speculum

S0074



Fishkind-Castroviejo Speculum

S0112



Lancaster Eye Speculum

CE

12

Speculum

S0120



Lieberman Speculum, V-wire

S0216.R




Lieberman Speculum, Reversible

S0118



Lieberman Speculum, K-wire

S0151



Schott Eye Speculum, Adjustable

S0189



Nasal Speculum

CE

13

Needle Holder

S0042



Castroviejo Needle Holder, Curved

S0046



Castroviejo Needle Holder Curved

S0043



Castroviejo Needle Holder with Lock

S0047



Castroviejo Needle Holder with Lock

S0044



Castroviejo Needle Holder Curved Delicate

S0045



Castroviejo Needle Holder Curved

S0018



Barraquer Needle Holder Curved Straight



Kalt Needle Holder



Callipers

S0140



Needle Holder
Curved
with lock

S0036



Castroviejo Caliper Curved

S0037



Castroviejo Caliper Straight

S0214



Helveston Scleral Ruler



Slit Blade



Implant Blade



Crescent Blade



MVR Blade



Lance Blade



Green Blade

Surgical Knives

Pk/10 per box
Sterile
Single use



15

Retractors

S0090



Helveston "Great Big Barbie"
Retractor, Concave

S0091



Helveston Barbie Retractor

S0092



Helveston Retractor, Concave
11mm

S0197



Freer Periosteal Elevator

S0213



S0212 Helveston Retractor,
7mm wide
S0213 Helveston Retractor,
9mm wide

S0109



S0109 Knapp Retractor, 4-prong, sharp
S0193 Stevenson Retractor 3x3 Prong
S0205 Knapp Retractor, 4-prong, blunt

S0055



S0055 Desmarres Lid Retractor 1
S0056 Desmarres Lid Retractor 2
S0057 Desmarres Lid Retractor 3

S0149



Schepens Orbital Retractor



Marker & Clamp

S0072




Fechtner Trabeculectomy Marker

S0050



Chayet Corneal Marker

S0080



Gass Scleral Marker

S0062



Downes Lid Clamp

S0122



Luntz-Dodick
Trabeculectomy Punch

S0152



Serrefine Clamp

S0152 Small
S0153 Large

S0103



Jones Towel Clamp

S0187



Backhaus Towel Clamp

CE

17

Hook

S0207



Helveston Finder Hook

S0208



Helveston Finder Hook, small

S0209



Helveston Muscle Hook, 10mm

S0210



Helveston Muscle Hook, 12mm

S0211



Helveston Muscle Hook, 8mm

S0215



Helveston Teaser Hook, 6mm

S0079



Gass Retinal Detachment Hook

S0075




Fixation Hook Double Sharp
Small



Probe & Rongeur

S0203



Lacrimal Probe

S0200



Pigtail Probe W/Suture Holes

S0086



Castroviejo Lacrimal Dilator

S0087



Harms Trabeculotomy Probe
Right
Left

S0190



Beyer Rongeur

CE S0188



Kerrison Rongeur

- S0188.1 Kerrison Rongeur, size 1
- S0188.1.5 Kerrison Rongeur, Size 1.5
- S0188.2 Kerrison Rongeur, size 2
- S0188.2.5 Kerrison Rongeur, Size 2.5
- S0188.3 Kerrison Rongeur, Size 3
- S0188.4 Kerrison Rongeur, Size 4

Chopper

S0049



Chang Phaco Chopper

S0105



Katzen LASIK Flap Unzipper

S0027



Blumenthal Conjunctiva Dissector

S0076



Flap Repositor

S0156



Sinskey II Lens Hook

S0104



Kansas Bisector

S0149




Schepens Orbital Retractor

S0157



Sinskey Lens Hook





- S0185 Micro-forcep (ILM, 23G)
- S0286 Micro-forcep (asymmetrical, 23G)
- S0287 Micro-forcep (asymmetrical, 25G)
- S0288 Micro-forcep (horizontal opening, serrated jaws, 23G)
- S0289 Micro-forcep (horizontal opening, end gripping jaws, 23G)
- S0294 Micro-forceps Horizontal Cutting (23G, 0.6mm)
- S0295 Micro-forceps Vertical cutting (23G, 0.6mm)
- S0291 Micro-forcep Intra Ocular Lens (IOL) handling angled
- S0159 Micro-forcep Serrated jaws forceps (23G, 0.6mm)
- S0159.S Micro-forcep Squeeze handle forceps -Titanium
- S0181 Micro Vitrectomy Scissor (Vertical Opening angle 60 degree 23G)



S0296 PVD Spatula (23G, 0.6mm)



BFN Flute Cannula Charles 23G 0.6mm





Alger Brush



Trepine

Cornea Instruments




Artificial Anterior Chamber




Endothelial Punch



Teflon Block




Suction Trepine



Vacuum Punch




S0155




**Simcoe (Irrigation-
aspiration Cannula)**
23G
Curved

S0195




Lacrimal Cannula
Malleable Tip- 23G
S0195 Curved
So195 Straight

S0096



Chu Hydrodissection cannula
25G
Flattened

S0276



Hydrodissection Cannula
27G
(Vertically flattened tip)

ACM-20



AC Maintainer
ACM-20 20G
ACM-23 23G
ACM-25 25G

CE 1118



Bimanual irrigation Hand Piece

Miscellaneous

S0156



Lens Manipulating Hook
Angled

S0182




Phaco Chopper Ball Tip

S0017



Iris Spatula
Round,
Smooth,
5mm)

S0077



Flieringa Rings

S0228



Ball Castray Point

S0242



Jaeger Lid Plate

S0051



Chuck Handle

S0227



Wright Facial Needle



Miscellaneous



Crash Cart
MI-39




LED Operation Theatre
Surgical Lights
994829



Mayo Surgical Trolley
A4S25



Vertical Autoclave
994489



Surgical
Instrument Tray
S0161



Sterilisation
Indicator Tape
916729



Eye Shields
Green 994908
Clear 916297
White 916298



Laser Eye Shields
LES-P-ND




Miscellaneous

S0184



Wells Enucleation Spoon

S0252



Vectis

S0252 Vectis Plain
S0253 Vectis Serrated

S0191



Mallet, Lead Filled

S0143




Paton Spatula & Spoon

S0199



Freer Lacrimal Chisel

S0198 Straight
S0199 Curved



Eye Drapes

60X60mm
90X110mm

CE

26

Forceps

Forcep tying straight with platform
Overall Length 87mm
ST0030



Forceps tying straight castroviejo with platform,
Over all Length 109mm
ST0030A



Forceps tying straight Kirby with platform, Over all Length 103mm
ST00333



Forceps tying curved with platform,
Over all Length 87mm
ST0036



Forceps tying curved castreivijo with platform,
Over all Length 109mm
ST0036A



Forcep tying curved Kirby with platform, Over all Length 103mm
ST0039



Forceps tying colibri with platform,
Over all Length 83mm
ST0042



Forceps tying colibri with platform,
Over all Length 83mm
ST0042



Forceps McPherson angled 10mm, with platform,
Over all Length 87mm
ST0048



Forceps McPherson Angled 12mm, castroviejo with platform,
Over all Length 109mm
ST0048A



Forceps McPherson angled 12mm, with platform,
Over all Length 85 mm
ST0054



Forceps

Forceps
capsulorhexis straight
overall length 85mm
ST0060



Utrata Forceps
capsulorhexis Castroviejo
overall length 109mm
ST0060A



Utrata Forceps capsulorhexis
kirby round handle
overall length 103mm
ST0063



Utrata Forceps
capsulorhexis curved shaft
kirby round handle
overall length 103mm
ST0063A



Utrata Forceps
capsulorhexis curved shaft
overall length 85mm
ST0066



Utrata Forceps
capsulorhexis curved shaft
castroviejo handle
overall length 109mm
ST0069



Capsulorhexis Forcep
23 gauge Curved Shaft
with micro grasping tip
ST0069A



Corneal Forceps toothed 0.12mm/0.2mm,
1x2 teeth with tying platform straight
overall length 87mm
ST0072



Corneal Forceps toothed
0.12mm/0.2mm, 1x2 teeth
with tying platform straight
castroviejo handle
overall length 107mm
ST0072A



Corneal Forceps toothed
0.12mm/0.2mm, 1x2 teeth
with tying platform kirby
round handle, overall length 103mm
ST0075



Corneal Forceps toothed 0.12/0.2mm,
1x2 teeth with tying platform curved
overall length 87mm
ST0078



Corneal Forceps toothed
0.12mm, 1x2 teeth with tying
platform/curved kirby round
handle, overall length 107mm
ST0078A



Forceps

Corneal Forceps colibri
0.12mm, 1x2 teeth with tying
platform kirby round handle
overall length 103mm
ST0081



Colibri Forceps 0.12mm,
1x2 teeth with tying platform
overall length 83mm
ST0084



Pollock Corneal Suturing forceps
very delicate double
1x2 teeth 0.12mm
overall length 83mm
ST0084A



Lims Forceps 0.12mm/0.2mm,
1x2 teeth with tying platform
overall length 93mm
ST0087



Lims Forceps pierse
type with tying platform
overall length 87mm
ST0087A



Lims Forceps 0.12mm/0.2mm,
1x2 teeth with tying platform
overall length 87mm
ST0087B



Forceps notched straight
pierce with tying platform
overall length 87mm
ST0090



Forceps notched curved
pierce with tying platform
overall length 87mm
ST0096



Forceps notched colibri
pierce with tying platform
overall length 83mm
ST00102



Forceps pierse type with
tying platform kirby round
handle, overall length 103mm
ST00105



Forceps capsule tag
1.5mm platform at tips
overall length 85mm
ST00108



Forceps

Forcep eye dressing
serrated delicate,
overall length 100mm
ST00111



Forcep fixation 1x2 teeth
overall length 90mm
ST00114



Forcep
superior rectus toothed
overall length 115mm
ST00117



Inomura
Capsulorhexis Forceps
overall length 115mm
ST00120



Chalozion
Hunt Forceps
overall length 95mm
ST00123



Chalozion Forceps
Oval Shape
overall length 95mm
ST00126



Jwellers Forceps
overall length 100mm
ST00129 Straight Fine
ST00132 Straight with platform



Faulkner Soft IOL
Holding Forceps
overall length 109mm
ST00135



Lens holding
forceps 'T' type
overall length 110mm
ST0003



Lens holding forceps
McPherson type
overall length 103mm
ST0009



Lens holding
forceps 'L' type
overall length 110mm
ST0015



30

Forceps/Needle Holder

Lens holding
forceps kratz
smooth curved jaws
overall length 110mm
ST0021



Lens holding
forceps Daljit model
overall length 85mm
ST0024



Needle holder barraquer
straight without lock
overall length 113mm
ST00300



Needle holder barraquer
straight with lock
overall length 113mm
ST00306



Needle holder barraquer
curved without lock,
overall length 113mm
ST00312



Needle Holder
Baraquer W/o Lock
Accurate Curved
overall length 113mm
ST00315



Needle holder barraquer
curved with lock,
overall length 113mm
ST00318



Needle holder
barraquer long model
overall length 135mm
ST00321



Needle holder
castroviejo curved
delicate without lock
overall length 113mm
ST00336



Needle Holder/Scissors/Speculum

Needle holder castroviejo
long curved
overall length 120mm
ST00342 with lock
ST00345 without lock



Needle holder castroviejo
long straight
overall length 120mm
ST00348 with lock
ST00351 without lock




Scissor corneal
universal curved
overall length 100mm
ST00240



Scissor corneal straight
overall length 100mm
ST00243



Scissor vannas curved
short blade/ long blade
overall length 95mm
ST00246



Scissor vannas straight
overall length 95mm
ST00249




Scissor vannas angled
upward capsulotomy
overall length 95mm
ST00252



Speculum, adjustable
overall length 82mm
ST00255 with 15mm solid blade for adults.
SST255A with 11mm solid blade for infants & child




Speculum adjustable spread
for phaco / lasik
overall size 82mm
ST00258




Lens Insert

Lieberman Speculum adjustable
open wire blade
ST00261



Speculum adjustable
open wire blade S.S.
ST00264

Kratz-Barraquer Wire Speculum
with open blades for unobstructed
access to the cataract incision adult / child
SS-267



Phaco silicon
soft lens folder
overall length 108mm
ST0021



Phaco silicon soft lens
insertor with lock
overall length 115mm
ST0027



Phaco acrylic lens folder
overall length 110mm
ST0030



Phaco acrylic lens insertor
overall length 108mm
ST0033



Phaco acrylic soft
lens folder box type
overall length 116mm
ST0036



Phaco acrylic soft lens
insertor to be used with
MP36 lens folder
overall length 112mm
ST0039



Folder Phaco acrylic soft IOL's
folding forceps with smooth platform
overall length 111mm
ST0042



Insertor for folding & holding the
shaft IOL's for insertion
overall length 100mm
ST0045



Injector/ Pre-chopper / Spatula / Repositor

Injector for acrylic soft lenses
overall length 121mm
ST0048



Pre-Chopper, Angled Round Handle
overall length 121mm
ST0054



Pre-Chopper, Curved Round Handle
overall length 121mm
ST0057



Akakoshi Pre-Chopper for cracking the soft or hard nucleus
overall length 121mm
ST0060



Repositor single ended malleable
overall length 123mm
ST00204



Spatula single ended
overall length 123mm
ST00207



Medelec repositor curved/angled, double ended,
overall length 145mm
ST00207A



Sinsky hook extra fine single ended
overall length 123mm
ST00210



Sinsky hook cum spatula double ended
overall length 145mm
ST00213



No hole IOL manipulator
overall length 115mm
ST00216



IOL manipulator ball ended
overall length 115mm
ST00219



Manipulator / Chopper / Spatula / Divider



Divider/Spatula/Irrigating Chopper/Y Rotator/Bimannual

Nucleus bisector for manual small incision surgery, to be used with MR255 overall length 110mm
ST00258



Nucleus cutting board for manual small incision surgery, to be used with MR255 & MR258, overall length 130mm
ST00261



Capsule polishing curette angled overall length 115mm
ST00264



Castroviejo Cyclodialysis Spatula Double ended Ultra Fine overall length 120mm
ST00267



Bi-manual Irrigating hand piece With S.S. 21G 2 irrigating ports overall length 115mm
ST00270



Bi-manual aspiration hand piece S.S. 21G with 0.3mm aspiration port. overall length 115mm
ST00273



Irrigating Nucleus Y-Shaped Manipulator 19G Shaft S.S. dual 0.6mm irrigating port Close to the tip for maximum irrigation. overall length 115mm
ST00276



Irrigating Nuclous Chopper with 0.6mm dual side port 20G overall length 115mm
ST00279 S.S. Gently Curved Shaft
ST00282 S.S. Angled Shaft



Fukasaku Hydro Chopper, blunt distal tip & front notch for engaging iris end. Opening 0.95mm O.D. overall length 115mm
ST00285



Nagahora Hydro Chopper Wedge shaped inferior edge with blunt distal round tip 1.5mm Long ended opening 0.95mm O.D. Overall length 115mm
ST00288



Rosen Hydro Chopper 1.1mm Long, Blunt end 0.95mm O.D. end opening. overall length 115mm
ST00291



Retractor / Vectis / Breaker

Desmarres Retractor
Lid Retractor, solid blade
Small / Med / Large
overall length 140mm
ST00294



Expressor cum
Vectis Double ended
overall length 150mm
ST00297



Expressor Cum
Spatula Double ended
overall length 150mm
ST00300



Lense Loop (Vectis)
overall length 140mm
ST00303



Wilder Lens Loop Serrated
overall length 140mm
ST00306



Catpow Retractor
overall length 130mm
ST00309



Blade breaker swiss made
overall length 115mm
ST00270



**Cataract Set with IOL Implantation Set
Set of 13 Instruments**



1. **ST0030** Tying forcep straight
2. **ST0036** Tying forcep curved
3. **ST0054** McPherson forcep angled
4. **ST0060** Capsulorhexis Forcep straight
5. **ST0087** Lims Forcep
6. **ST00114** Fixation Forcep
7. **ST00117** Rectus Forcep toothed superior
8. **ST00312** Barraquer Needle Holder
9. **ST00240** Corneal Universal Scissor Curved
10. **ST00204** Repositor Single ended
11. **ST00207A** Repositor cum Spatula double ended
12. **ST00210** Sinsky Hook extra fine single ended
13. **ST00219** I.O.L. Manipulator ball ended

Micro Surgery Sets

**Cataract Set
Set of 6 Instruments**



- | | |
|---------------|----------------|
| ST0030 | ST00312 |
| ST0054 | ST00240 |
| ST0087 | ST00210 |

**Cataract Set
Set of 4 Instruments**



- | | |
|---------------|---------------|
| ST0030 | ST0087 |
| ST0054 | SST312 |

**Phaco Set
Set of 4 Instruments**



- | | |
|---------------|----------------|
| ST0030 | ST00222 |
| ST0033 | ST00243 |



Bandage Contact Lens



The aims of a bandage contact lens are to:

- Relieve some causes of eye pain
- Promote healing of the cornea
- Mechanical protection and support
- Maintenance of moist corneal surface

The can also help in stabilizing wounds, cuts and holes in the cornea by reducing pain and preventing secondary corneal infection by acting as a barrier to microbes (bacteria/germs).

*Customized Bandage Contact Lenses are also available

- Extended Wear
- High Water Content
- Optimum Oxygen permeability

A bandage contact lens protects an injured or diseased cornea from the mechanical rubbing of the blinking eyelids, therefore allowing it to heal. It is a normal (usually soft, but not always) contact lens has no power. It is not intended to improve your vision.

Code	Products	Material	Sterile	Packing
BCL142/20	BCL 14.2	HEMA	EO	Pack/20
BCL160/20	BCL 16mm	HEMA	EO	Pack/20
BCL18/20	BCL 18mm	HEMA	EO	Pack/20
BCL20/20	BCL 20mm	HEMA	EO	Pack/20
BCL142/6	BCL 14.2	HEMA	EO	Pack/6
BCL160/6	BCL 16mm	HEMA	EO	Pack/6
BCL18/6	BCL 18mm	HEMA	EO	Pack/6
BCL20/6	BCL 20mm	HEMA	EO	Pack/6

Eye Sphere



When an eye is severely damaged by trauma or disease, it may be necessary to remove it surgically. The most common procedures are evisceration (removal of the contents of the eye) or enucleation (removal of the eye itself).

Eye spheres are permanent implants to occupy the cavity that results from surgery. Conformers are temporary devices, which maintain the shape of the eye and prevent closure or adhesion during the healing process. When healing is complete, conformers are replaced with a prosthesis.


AOPL0020/U/25/10	Eye Sphere	PMMA	10mm	Unsterile	Pack/25
AOPL0020/U/25/12	Eye Sphere	PMMA	12mm	Unsterile	Pack/25
AOPL0020/U/25/14	Eye Sphere	PMMA	14mm	Unsterile	Pack/25
AOPL0020/U/25/16	Eye Sphere	PMMA	16mm	Unsterile	Pack/25
AOPL0020/U/25/18	Eye Sphere	PMMA	18mm	Unsterile	Pack/25
AOPL0020/U/25/20	Eye Sphere	PMMA	20mm	Unsterile	Pack/25
AOPL0020/U/25/22	Eye Sphere	PMMA	22mm	Unsterile	Pack/25
AOPL0020/S/25/24	Eye Sphere	PMMA	24mm	Unsterile	Pack/25
AOPL0020/S/10	Eye Sphere	PMMA	10mm	Sterile	Pack/1
AOPL0020/S/12	Eye Sphere	PMMA	12mm	Sterile	Pack/1
AOPL0020/S/14	Eye Sphere	PMMA	14mm	Sterile	Pack/1
AOPL0020/S/16	Eye Sphere	PMMA	16mm	Sterile	Pack/1
AOPL0020/S/18	Eye Sphere	PMMA	18mm	Sterile	Pack/1
AOPL0020/S/20	Eye Sphere	PMMA	20mm	Sterile	Pack/1
AOPL0020/S/22	Eye Sphere	PMMA	22mm	Sterile	Pack/1
AOPL0020/S/24	Eye Sphere	PMMA	24mm	Sterile	Pack/1
AOPL0020/SS/12	Eye Sphere	Silicone	12mm	Sterile	Pack/1
AOPL0020/SS/14	Eye Sphere	Silicone	14mm	Sterile	Pack/1
AOPL0020/SS/16	Eye Sphere	Silicone	16mm	Sterile	Pack/1
AOPL0020/SS/18	Eye Sphere	Silicone	18mm	Sterile	Pack/1
AOPL0020/SS/20	Eye Sphere	Silicone	20mm	Sterile	Pack/1

Carter Eye Sphere Introducer



When an eye is severely damaged by trauma or disease, it may be necessary to remove it surgically. The most common procedures are evisceration (removal of the contents of the eye) or enucleation (removal of the eye itself).

Eye spheres are permanent implants to occupy the cavity that results from surgery. Conformers are temporary devices, which maintain the shape of the eye and prevent closure or adhesion during the healing process. When healing is complete, conformers are replaced with a prosthesis.



Orbital Implant Sizer



Product Code: AOPL0021

Akriti sizers are designed to accurately select the appropriate orbital implant diameter in case of enucleation / evisceration.

Main characteristics:

- 5 sizes: 12,14, 16, 18 and 20 mm
- Made of PMMA
- Sterile, Single use

AOPL0021-12	Implant sizer 12mm	Clear	20 mm	PMMA	Box of 1
AOPL0021-14	Implant sizer 14mm	Clear	23 mm	PMMA	Box of 1
AOPL0021-16	Implant sizer 16mm	Clear	26 mm	PMMA	Box of 1
AOPL0021-18	Implant sizer 18mm	Clear	20 mm	PMMA	Box of 1
AOPL0021-20	Implant sizer 20mm	Clear	23 mm	PMMA	Box of 1

Custom Packing Available with MOQ

Symblepharon Ring

Symblepharon rings are designed to prevent retraction of the conjunctival cul-de-sacs without touching the cornea.



Main characteristics:

- Highly polished and atraumatic edges
- Translucent/Clear
- Available in 3 sizes
- Sterile, Single use
- No More Steep Curve open fornix conformers , which
- make your corneal graft difficult to stay
- Precisely made
- **Size**
- Small. : 20mm
- Medium : 23mm
- Large : 25mm
- Sterile



AOPL0023/1	Symblepharon Ring	Small	20 mm	PMMA	Box of 1
AOPL0023/2	Symblepharon Ring	Medium	22 mm	PMMA	Box of 1
AOPL0023/3	Symblepharon Ring	Large	25 mm	PMMA	Box of 1
AOPL0023/4	Symblepharon Ring	Small	20 mm	Polyamide	Box of 1
AOPL0023/5	Symblepharon Ring	Medium	22 mm	Polyamide	Box of 1
AOPL0023/6	Symblepharon Ring	Large	25 mm	Polyamide	Box of 1

Custom Packing Available with MOQ



Ocular Conformer

Perforated conformers from Akriti maintain the conjunctival cul-de-sacs after enucleation or evisceration without hindering lid closure.

Main characteristics:

- Perforations allow for postoperative instillation of eye drops and drainage of secretions
- Available in 2 models: rigid or flexible
- Translucent
- Sterile, Single use
- No More Steep/ Flat Curve for your ocular conformers , which make your tissue aaposition opens and make it implant rejection



Size:

- Small 20 mm
- Medium 23mm
- Large 26mm
- Sterile
- CE

AOPL0022-001	Perforated Conformers	Rigid & small	20 mm	PMMA	Box of 1
AOPL0022-002	Perforated Conformers	Rigid & medium	23 mm	PMMA	Box of 1
AOPL0022-003	Perforated Conformers	Rigid & large	25 mm	PMMA	Box of 1
AOPL0022-004	Perforated Conformers	Flexible & small	20 mm	Polyamide 12	Box of 1
AOPL0022-005	Perforated Conformers	Flexible& medium	23 mm	Polyamide 12	Box of 1
AOPL0022-006	Perforated Conformers	Flexible & large	25 mm	Polyamide 12	Box of 1

Custom Packing Available with MOQ



Artificial Eyes


The **Artificial Eye** is fitted to replace the lost volume of the socket and to restore natural appearance and movement. Artificial eye is fitted after enucleation of the anterior segment of ocular evisceration. It is made of acrylic and it can withstand time, climate or any other element that could affect the eye. Artificial eye also readily adapts itself to the temperature of the socket.

The object of fitting an ocular prosthesis (artificial eye) is:

- To maintain shape and the movement of the eyelid and for cosmetic purposes
- To retain the tone of facial muscles to prevent alterations of features
- To prevent foreign bodies from entering the ocular cavity
- To prevent accumulation of fluid in the cavity and to direct secretions to proper channels.





















Available Sizes : Right and Left Eye; 20mm, 22mm, 24mm, 26mm and 28mm.

Pupil	Black
Material	Acrylic Powder
Size	15MM TO 28 MM
Color	Black, Brown, Light Blue, Dark Blue, Green



Product Code: AE

Click on the images below to select your Eye

Sclera	Iris						Image
 LOC-LV	Light Brown			Blue			 Iris Id: <input type="text" value="lbs3"/> Sclera Id: <input type="text" value="mlw-mv"/>
							
	LBS3	LBS5	LBS7	BLS1	BLS4	BLS5	
	Medium Brown			Grey			
 MLW-MV							
	MBS1	MBS4	MBS6	GYS1	GYS2	GYS3	
 C-DV	Dark Brown			Green			
							
	DBS2	DBS4	DBS7	GNS1	GNS3	GNS4	
	Black						
							
	DBS3	BKS2					

Corneal Shields



Main characteristics:

- To protect the eye during lacrimal or lid surgery
- Supplied with a latex suction cup for easy placement & removal
- Not for laser treatment
- Sterile, Single use

AOPL0024/1	Ocular shield	Small	20 mm	PMMA	Box of 1
AOPL0024/2	Ocular shield	Medium	22 mm	PMMA	Box of 1
AOPL0024/3	Ocular shield	Large	25 mm	PMMA	Box of 1

Suction Cups



Artificial **Eye Suction Cups** aid in the removal and placement of artificial or **prosthetic eyes**. They have thin walls for easier squeezing and are washable with soap and water.

Corneal Shields



Main characteristics:

- To protect the eye during lacrimal or lid surgery
- Supplied with a latex suction cup for easy placement & removal
- Not for laser treatment
- Sterile, Single use

AOPL0024/1	Ocular shield	Small	20 mm	PMMA	Box of 1
AOPL0024/2	Ocular shield	Medium	22 mm	PMMA	Box of 1
AOPL0024/3	Ocular shield	Large	25 mm	PMMA	Box of 1

Gold Weight Implant



Implantable eyelid gold weights are used for the treatment of lagophthalmos. Thinprofile eyelid implants are 40 percent thinner than the original Contour eyelid implants. This provides a sleek, spherical curvature with tapered edges and smooth, rounded corners. For added freedom in suturing position, channel-free 1.0 mm holes are used. Although the time-tested AKRITI eyelid implants are still a highly effective treatment option, ThinProfile implants offer a more aesthetically pleasing treatment for lagophthalmos.


Code	Products	Material	Thickness	Sterile	Packing
AOPL006	Eye Lid Gold weight 0.6gm	24K Gold	1.0 mm	EO	Box of 1
AOPL008	Eye Lid Gold weight 0.8gm	24K Gold	1.0 mm	EO	Box of 1
AOPL010	Eye Lid Gold weight 1.0	24K Gold	1.0 mm	EO	Box of 1
AOPL012	Eye Lid Gold weight 1.2gm	24K Gold	1.0 mm	EO	Box of 1
AOPL014	Eye Lid Gold weight 1.4gm	24K Gold	1.0 mm	EO	Box of 1
AOPL016	Eye Lid Gold weight 1.6gm	24K Gold	1.0 mm	EO	Box of 1
AOPL018	Eye Lid Gold weight 1.8gm	24K Gold	1.0 mm	EO	Box of 1
AOPL020	Eye Lid Gold weight 2.0gm	24K Gold	1.0 mm	EO	Box of 1
AOPL022	Eye Lid Gold weight 2.2gm	24K Gold	1.0 mm	EO	Box of 1
AOPL024	Eye Lid Gold weight 2.4gm	24K Gold	1.0 mm	EO	Box of 1
AOPL026	Eye Lid Gold weight 2.6gm	24K Gold	1.0 mm	EO	Box of 1
AOPL028	Eye Lid Gold weight 2.8gm	24K Gold	1.0 mm	EO	Box of 1

Lacrimal Intubation Sets

Lacrimal Intubation set comprising of stainless steel probes with silicon tubing length 35 cm. Packed in ready to use sterile packs.


Sizes:

- 11 cm x 23 G
- 7.5 cm x 23 G
- 11 cm x 23 G with Olive tips
- 7.5 cm x 23 G with Olive tips
- 11 cm x 27 G



Code	Products	Material	Ste rile	Packing
LIS23G7C	Lacrimal Intubation Set 23G 7.5 cm	Silicone	EO	Pack/10
LIS23G11C	Lacrimal Intubation Set 23G 11 cm	Silicone	EO	Pack/10
LICO27G11	Lacrimal Intubation Set Olive tip 27G 11 cm	Silicone	EO	Pack/10
LICO27G11	Lacrimal Intubation Set Olive tip 23G 11 cm	Silicone	EO	Pack/10
RHD 1	Retrival Hook Devices		EO	Pack/10

Lester Johns Tube With Flange



Dacryocystorhinostomy (DCR) surgery provides relief to millions of people every year who suffer from tear duct problems ranging from collapsed tear ducts and excessive tearing caused by accident or injury, birth defects or environmental strain. Left untreated, patients may experience swelling, infections, increased swelling, tearing ("wet eye") and severe pain.


In 1961, Lester T. Jones and Gunther Weiss paired up to develop a solution for tear duct problems. They developed highly specialized glass tubes (called "Lester T. Jones Tear Duct Tubes") designed for insertion into the nasolacrimal duct. The body develops scar tissue around the tube, which holds it in place, forming a new permanent tear duct.

3.0mm

Code	Products	Material	Ste rile	Packing
1000-3510	Lester Johns Tube3.5X10 mm With Flange	Pyrex Glass	No	Pack/1
1000-3511	Lester Johns Tube3.5X11 mm With Flange	Pyrex Glass	No	Pack/1
1000-3512	Lester Johns Tube3.5X12 mm With Flange	Pyrex Glass	No	Pack/1
1000-3513	Lester Johns Tube3.5X13 mm With Flange	Pyrex Glass	No	Pack/1
1000-3514	Lester Johns Tube3.5X14 mm With Flange	Pyrex Glass	No	Pack/1
1000-3515	Lester Johns Tube3.5X15 mm With Flange	Pyrex Glass	No	Pack/1
1000-3516	Lester Johns Tube3.5X16 mm With Flange	Pyrex Glass	No	Pack/1
1000-3517	Lester Johns Tube3.5X17 mm With Flange	Pyrex Glass	No	Pack/1
1000-3518	Lester Johns Tube3.5X18 mm With Flange	Pyrex Glass	No	Pack/1
1000-3520	Lester Johns Tube3.5X20 mm With Flange	Pyrex Glass	No	Pack/1

50

Lester Johns Tube With Flange



Dacryocystorhinostomy (DCR) surgery provides relief to millions of people every year who suffer from tear duct problems ranging from collapsed tear ducts and excessive tearing caused by accident or injury, birth defects or environmental strain. Left untreated, patients may experience swelling, infections, increased swelling, tearing ("wet eye") and severe pain.


In 1961, Lester T. Jones and Gunther Weiss paired up to develop a solution for tear duct problems. They developed highly specialized glass tubes (called "Lester T. Jones Tear Duct Tubes") designed for insertion into the nasolacrimal duct. The body develops scar tissue around the tube, which holds it in place, forming a new permanent tear duct.

3.5mm

Code	Products	Material	Ste rile	Packing
1000-3510	Lester Johns Tube3.5X10 mm With Flange	Pyrex Glass	No	Pack/1
1000-3511	Lester Johns Tube3.5X11 mm With Flange	Pyrex Glass	No	Pack/1
1000-3512	Lester Johns Tube3.5X12 mm With Flange	Pyrex Glass	No	Pack/1
1000-3513	Lester Johns Tube3.5X13 mm With Flange	Pyrex Glass	No	Pack/1
1000-3514	Lester Johns Tube3.5X14 mm With Flange	Pyrex Glass	No	Pack/1
1000-3515	Lester Johns Tube3.5X15 mm With Flange	Pyrex Glass	No	Pack/1
1000-3516	Lester Johns Tube3.5X16 mm With Flange	Pyrex Glass	No	Pack/1
1000-3517	Lester Johns Tube3.5X17 mm With Flange	Pyrex Glass	No	Pack/1
1000-3518	Lester Johns Tube3.5X18 mm With Flange	Pyrex Glass	No	Pack/1
1000-3520	Lester Johns Tube3.5X20 mm With Flange	Pyrex Glass	No	Pack/1

51

Lester Johns Tube With Flange



Dacryocystorhinostomy (DCR) surgery provides relief to millions of people every year who suffer from tear duct problems ranging from collapsed tear ducts and excessive tearing caused by accident or injury, birth defects or environmental strain. Left untreated, patients may experience swelling, infections, increased swelling, tearing ("wet eye") and severe pain.

In 1961, Lester T. Jones and Gunther Weiss paired up to develop a solution for tear duct problems. They developed highly specialized glass tubes (called "Lester T. Jones Tear Duct Tubes") designed for insertion into the nasolacrimal duct. The body develops scar tissue around the tube, which holds it in place, forming a new permanent tear duct.

4.0mm

Code	Products	Material	Ste rile	Packing
1000-4010	Lister Johns Tube 4.0X 10mm	Pyrex Glass	No	Pack/1
1000-4011	Lister Johns Tube 4.0X 11mm	Pyrex Glass	No	Pack/1
1000-4012	Lister Johns Tube 4.0X 12mm	Pyrex Glass	No	Pack/1
1000-4013	Lister Johns Tube 4.0X 13mm	Pyrex Glass	No	Pack/1
1000-4014	Lister Johns Tube 4.0X 14mm	Pyrex Glass	No	Pack/1
1000-4015	Lister Johns Tube 4.0X 15mm	Pyrex Glass	No	Pack/1
1000-4016	Lister Johns Tube 4.0X 16mm	Pyrex Glass	No	Pack/1
1000-4017	Lister Johns Tube 4.0X 17mm	Pyrex Glass	No	Pack/1
1000-4018	Lister Johns Tube 4.0X 18	Pyrex Glass	No	Pack/1
1000-4020	Lister Johns Tube3.5 20	Pyrex Glass	No	Pack/1

Universal Eye Shields



Universal, clear, lightweight plastic with smooth edges & ventilation holes for patient comfort.

- Universal design fits both left and right eyes, thus minimising inventory
- Transparent material allows limited vision if required
- High quality smooth, ergonomic design for optimum patient comfort


AOPL0025/U/50	Ocular Eye Shield	A	Unsterile	Polycarbonate	Box of 100
AOPL0025/U/500	Ocular Eye Shield	A	Unsterile	Polycarbonate	Box of 500
AOPL0025/S/50	Ocular Eye Shield	A	Sterile	Polycarbonate	Box of 50
AOPL0028/U/100	Ocular Eye Shield	C	Unsterile	Polycarbonate	Box of 100
AOPL0028/U/1000	Ocular Eye Shield	C	Unsterile	Polycarbonate	Box of 500
AOPL0028/S/50	Ocular Eye Shield	C	Sterile	Polycarbonate	Box of 50



i-SHIELD

Custom Packing Available with MOQ

Eye Shields - Green



Eye shields are used after surgery to protect the eye and hold dressings in place. They are made of clear shatterproof plastic in a universal shape for use over either eye. Holes in the eye shield provide ventilation.

AOPL0026/U/100/G	Ocular Eye Shield	A	Unsterile	Polypropylene	Box of 100
AOPL0026/U/500/G	Ocular Eye Shield	A	Unsterile	Polypropylene	Box of 500
AOPL0026/S/50/G	Ocular Eye Shield	A	Sterile	Polypropylene	Box of 50
AOPL0027/U/100	Ocular Eye Shield	C	Unsterile	Polypropylene	Box of 100
AOPL0027/U/1000	Ocular Eye Shield	C	Unsterile	Polypropylene	Box of 500
AOPL0027/S/50	Ocular Eye Shield	C	Sterile	Polypropylene	Box of 50



Colors Available



G

W


B

Y

i-SHIELD

Custom Packing Available with MOQ

i-SAFE Clinical Safety Glasses




- Breathable
- Optical Clear
- Shatter proof
- Can be worn over the Glasses
- Perfect for fitting over most prescription lens
- Added feature of vents that keep the glasses from fogging
- The goggles are uv protected.impact resistant lens, side-shield protection
- Can be used while riding motorcycles, Impact Resistant Lens, Also Be Used As Bike Driving Goggles, Crystal Clear Vision

BS EN166: 2002 Compliance


AOPL0014CL	Clinical Safety	I-SAFE Medical Safety Glasses	Polycarbonate	Box of 20
AOPL0014YL	Night Vision	I-SAFE Night Vision Glasses	Polycarbonate	Box of 20
AOPL0014 BK	Post operative Glasses	I-SAFE Medical Safety Glasses Post Operative	Polycarbonate	Box of 20

55


Circling Bands



To be used with 1, 2, 4.




To be used with 1, 3, 4.




To be used with 1, 3, 4.

Silicone Strips



To be used with 5, 6, 7.







To be used with 5, 6, 7.

Silicone Tires




Concave



Convex



Asymmetrical



Retina Tire and Retina Band made out of Implantable Medical Grade Silicon Rubber. This is a temporary implant used to close the retinal holes or to reattach the detached Retina and relieve vitreous traction.

The Retina Tire is available mainly in three different styles. Retina Tire Style 279 and 281 are symmetrical while Retina Tire Style 280 is asymmetrical.

The Retina Band is available in Style 240, 241

Retinal Sleeve available for 240, 270, 287

RB240	Retinal Band 240	Silicone	125X2.5X0.6	Symmetrical	Box of 10
RB241	Retinal Band 241	Silicone		Symmetrical	Box of 10
RT279	Retinal Tire 279	Silicone	2.5X 5.0mm	Symmetrical	Box of 10
RT280	Retinal Tire 280	Silicone	2.5X 7.0mm	Symmetrical	Box of 10
RT281	Retinal Tire 281	Silicone	2.5X 7.0mm	Asymmetrical	Box of 10
RS240	Retinal Band Sleeve 240	Silicone		Symmetrical	Box of 10
RS270	Retinal Band Sleeve 270	Silicone		Symmetrical	Box of 10
RS287	Retinal Band Sleeve 287	Silicone		Symmetrical	Box of 10

Eye Drapes



Ophthalmic Drapes are used extensively for covering the body parts except the eye area during operation.

Disposable Ophthalmic Drape with Pouch

Item Code	Size
001002-50x40	40 cm x 50 cm
001002-90X75	90 cm x 75 cm
001002-80X100	80 cm x 100 cm
001002-100x100	100 cm x 100 cm
001002-100X120	100 cm x 120 cm

MOQ	5000 Piece
Usage/Application	Hospital
Brand	Akriti
Material	Poly Material
Sterility	Sterile
Quantity Per Pack	50
Certification	ISO 13485-2016
Fluid Impervious	Yes

India's 1st Microprocessor Controlled Radio Frequency Surgical Unit

SKU 0014



Features :


- ❖ Output Frequency 2.0 Mhz
- ❖ Line Voltage 110V/220V/240V
- ❖ Line Frequency 50-60Hz

Output Power :

- ❖ RMS 90 Watts
- ❖ Peak 120 Watts
- ❖ Weight 1.2 Kg

Benefits to your practice :

- ❖ User friendly design very easy to operate
- ❖ Cost Effective, disposable hand pieces
- ❖ Foot switch activated with Finger switch control
- ❖ 50-60 Microns thermal heat spread



Set of blades to be used

Features :

- ❖ Precise Incisions
- ❖ Versatile
- ❖ Quick recovery for tissue reconstruction
- ❖ Low post operative pain to Patient
- ❖ Minimum Burning of tissue
- ❖ Min Smoke and plume, which allows better Visualization while reducing the odor



Sterile disposable low cost hand pieces

Compliance: **ISO 9001 :2008**
ISO 9001 :13485

SKU 0013



Portable Electrolysis Unit

- ❖ Line Voltage 110V/220V/240V
- ❖ Line Frequency 50-60Hz

MANUAL PHACOFRACTURE CATARACT EXTRACTION



AOLSSI- MANUAL PHACOFRACTURE CATARACT EXTRACTION

as used by Eduardo Alfonso, MD

- | | | | |
|-------------------|------------------------------------|---------------------|---|
| 1. S0120 | Liebe rman Specu l u m, V-wire | 15. S0032 | Bucher-Kansas Nucleus Trisector |
| 2. S0069 | Eye Scissors, standard, straight | 16. S0104 | Kansas Nucleus Bisector |
| 3. S0186 | Westcott Tenotomy Scissors | 17. S0008 | Alfonso Forceps, 3x4 teeth |
| 4. S0040 | Castroviejo Forceps, 0.12mm | 18. S0012 | Alfonso-McIntyre Irrigating Spoon, sma |
| 5. S0051 | Chuck Hand e (3) | 19. S0011 | Alfonso-McIntyre Irrigating Spoon, large |
| 6. S0146 | Performance Lamelar Blade #66 (bx) | 20. S0127 | McIntyre Infusion Handle |
| 7. MM-1500 | Performance Super Sharp Blade (bx) | 21. MM051-27 | Air njection Cannula, 27ga. |
| 8. MM-5570 | Performance Phaco Blade (bx) | 22. S0013 | Anis Cortex Aspirating Cannula, straight |
| 9. S0037 | Castroviejo Caliper, straight | 23. S0031 | Brierley Chamber Maintainer |
| 10. S0177 | Utrata Capsulorrhexis Forceps | 24. S0024 | Bipolar Forceps, straight, pointed |
| 11. S0111 | Kuglen Lens Manipulator, ang ed | 25. S0061 | Dodick IOL Folding Forceps |
| 12. S0006 | Alfonso Cutting Platform, 2mm | 26. S0035 | Buratto Soft IOL Inserter |
| 13. S0007 | Alfonso Cutting Platform, 3mm | 27. S0154 | Silicone Tubing w/adaptor for K7-6750 |
| 14. S0009 | Alfonso Nucleus Trisector, oval | 28. S0161 | Sterilizing Case, plastic (not illustrated) |

* not illustrated

Technique Video Available


PLANNED EXTRA-CAPSULAR CATARACT EXTRACTION



AOLSS2- PLANNED EXTRA-CAPSULAR CATARACT EXTRACTION

- | | | | | | |
|----|----------------|---|----|-----------------|--|
| 1 | S0025 | Bishop-Harmon A/C Irrigator | 23 | S0133 | McPherson Tying Forceps, straight |
| 2 | S0103 | Jones Towel Clamp | 24 | S0132 | McPherson Tying Forceps, angled |
| 3 | S0089 | Hartmann Mosquito Forceps, straight | 25 | S0106 | Kelman-McPherson Forceps |
| 4 | S0152 | Serrefine, small, straight | 26 | S0177 | Utrata Capsulorrhexis Forceps |
| 5 | S0118 | Lieberman Speculum, K-wire | 27 | S0156 | Sinsky II Lens Hook, angled |
| 6 | S0019 | Barraquer Wire Speculum, large | 28 | S0110 | Kuglen Iris Hook, angled |
| 7 | S0051 | Chuck Handle (2) | 29 | S0017 | Barraquer Iris Spatula |
| 8 | MM-1500 | Performance Super Sharp Blade, 15 deg., 5.0mm | 30 | S0018 | Barraquer Needle Holder, curved |
| 9 | S0144 | Performance Lamellar Blade #66 | 31 | S0046 | Castroviejo Needle Holder, straight |
| 10 | S0059 | Diamond Step Knife | 32 | S0116 | Lewis Lens Loop |
| 11 | PV-S | PVA Sponges box of 50 | 33 | S0010 | Alfonso Trisector, oval |
| 12 | S0039 | Castroviejo Corneal Scissors, small, right | 34 | S0005 | Alfonso Cutting Platform |
| 13 | S0038 | Castroviejo Corneal Scissors, small, left | 35 | S0008 | Alfonso Forceps, 3x4 teeth |
| 14 | S0186 | Westcott Tenotomy Scissors | 36 | S0012 | Alfonso-McIntyre Irrigating Spoon, small |
| 15 | S0134 | Westcott Stitch Scissors | 37 | S0011 | Alfonso-McIntyre Irrigating Spoon, large |
| 16 | S0179 | Vannas Scissors, curved | 38 | S0014 | Anis Cortex Aspiration Cannula |
| 17 | S0163 | Stern-Gills Scissors | 39 | MM051-27 | Air Injection Cannula, 27ga. |
| 18 | S0036 | Castroviejo Caliper | 40 | MM051-26 | McIntyre A/C Cannula, 26ga. |
| 19 | S0054 | Conjunctiva Forceps | 41 | S0101 | Jensen Polisher, 23ga. |
| 20 | S0088 | Harms-Colibri Forceps | 42 | S0155 | Simcoe I/A Cannula, curved |
| 21 | S0029 | Bonn Forceps, 0.12mm | 43 | S0031 | Brierley Chamber Maintainer |
| 22 | S0026 | Bishop-Harmon Forceps | 44 | S0161 | Sterilizing Case, double (not shown) |

LACRIMAL SURGERY



AOLSS3- LACRIMAL SURGERY

1	S0025	Bishop-Harmon A/C Irrigator	17	S0199	Freer Lacrimal Chisel, curved
2	S0187	Backhaus Towel Clamp (2)	18	S0200	Pigtail Probe w/Suture Holes
3	S0085	Halsted Mosquito Forceps, straight (2)	19	S0201	Crawford Retrieving Hook
4	S0084	Halsted Mosquito Forceps, curved (2) (not shown)	20	S0202	Crawford Lacrimal Intubation Set
5	S0188	Kerrison Rongeur, size 1	21	S0203	Bowman Lacrimal Probe Set
6	S0189	Nasal Speculum, adult	22	S0063	Dressing Forceps, delicate, straight
7	S0190	Beyer Rongeur	23	S0204	Sauer Fixation Forceps, 1x2 teeth
8	S0191	Mallet, lead filled stainless steel	24	S0028	Bonaccolto Conjunctiva Forceps
9	S0192	Performance Miniature Blade #65	25	S0041	Castroviejo Forceps, 0.3mm
10	S0051	Chuck Handle	26	S0205	Knapp Retractor, 4-prong, blunt (2)
11	S0193	Stevenson Retractor 3x3 Prong	27	S0164	Stevens Scissors, curved
12	S0194	Lacrimal Cannula, straight	28	S0070	Eye Scissors, straight
13	S0195	Lacrimal Cannula, curved	29	S0043	Castroviejo Needle Holder, curved w/lock
14	S0196	Castroviejo Lacrimal Dilator	30	S0134	Westcott Stitch Scissors
15	S0197	Freer Periosteal Elevator	31	PV-S	PVA Sponges box of 50 (not illustrated)
16	S0198	West Lacrimal Chisel, straight	32	S0161	Sterilizing Case, plastic (not illustrated)


PLANNED EXTRA-CAPSULAR CATARACT EXTRACTION



AOLSS4- EXTRA-CAPSULAR CATARACT EXTRACTION USING PHACOEMULSIFICATION

- | | | | |
|------------|-------------------------------------|----------------|---|
| 1. S0025 | Bishop-Harmon A/C Irrigator | 22. S0177 | Ultracapsulorrhexis Forcep |
| 2. S0118 | Lieberman Speculum, K-wire | 23. S0107 | Kelman-McPherson Forceps, angled |
| 3. S0167 | Temporal Speculum, open wire | 24. S0066 | Ernest-McDonald IOL Inserter |
| 4. S0058 | Diamond Knife, tri-facet | 25. S0121 | Livernois IOL Folder |
| 5. K2-6607 | Phaco-4 Diamond Knife 3.0mm, angled | 26. S0140 | Katena Needle Holder, curved, w/loc |
| 6. S0051 | Chuck Hand e | 27. S0065 | Ernest Nucleus Cracker |
| 7. MM-1500 | Performance Super Sharp Blade (bx) | 28. S0108 | Keuch Pupil Dilator |
| 8. S0089 | Hartmann Mosquito Forceps, straight | 29. S0155 | Simcoe I/A Cannula, curved |
| 9. S0103 | Jones Towel Clamp | 30. S0130 | McIntyre-Binkhorst Cannula, right |
| 10. S0148 | Rosen Phaco Splitter | 31. S0129 | McIntyre-Binkhorst Cannula, left |
| 11. S0110 | Kuglen ris Hook, angled | 32. MM-1073-25 | Hydrodissection Cannula, 25ga. |
| 12. S0100 | Jaffe-Bechert Rotator, 9mm | 33. MM051-27 | Air Injection Cannula, 27ga. |
| 13. S0114 | Lester IOL Manipulator, angled | 34. S0101 | Jensen Capsule Polisher, 23ga. |
| 14. S0186 | Westcott Tenotomy Scissors | 35. S0166 | Stevens Sub-Tenon's Cannula, 19ga. |
| 15. S0037 | Castroviejo Caliper, straight | 36. S0073 | Fine-Thornton Fixation Ring, 13mm |
| 16. S0162 | Stern-Gills Capsulotomy Scissors | 37. S0023 | Bipoar Forceps, angled, blunt |
| 17. S0053 | Colibri Forceps, 0.12mm teeth | 38. S0022 | Bipoar Cord |
| 18. S0040 | Castroviejo Forceps, 0.12mm | 39. S0128 | McIntyre Reverse Cystotome |
| 19. S0071 | Fechtner Conjunctiva Forceps | 40. S0179 | Vannas Scissors, curved |
| 20. S0169 | Tennant Tying Forceps, straight | 41. PV-S | PVA-Sponge Spears, (bx) (not illustrated) |
| 21. S0168 | Tennant Tying Forceps, curved | 42. S0161 | Sterilizing Case, plastic (not illustrated) |


VITREORETINAL SURGERY



AOLSS5 - VITREORETINAL SURGERY

1	S0074	Fishkind-Castroviejo Speculum	19	S0183	Watzke Sleeve Spreading Forceps
2	S0152	Serrefine, small, straight (4)	20	S0028	Bonaccolto Conjunctiva Forceps
3	S0085	Halsted Mosquito Forceps, straight (2)	21	S0141	Nugent Utility Forceps, serrated
4	S0070	Eye Scissors, straight	22	S0040	Castroviejo Forceps, 0.12mm
5	S0025	Bishop-Harmon A/C Irrigator	23	S0041	Castroviejo Forceps, 0.3mm
6	S0165	Stevens Scissors, straight	24	S0026	Bishop-Harmon Forceps (2)
7	S0164	Stevens Scissors, curved	25	S0043	Castroviejo Needle Holder, curved w/lock
8	S0186	Westcott Tenotomy Scissors	26	S0047	Castroviejo Needle Holder, straight w/lock
9	S0134	Westcott Stitch Scissors	27	S0133	McPherson Tying Forceps, straight
10	S0145	Performance Miniature Blade #64 (bx)	28	S0091	Helveston Barbie Retractor
11	S0051	Chuck Handle	29	S0092	Helveston Retractor, concave, 11mm
12	MM-5563	Performance MVR Knife, straight, 19ga.	30	S0180	Viscoelastic Aspiration Cannula, 22ga.
13	S0149	Schepens Orbital Retractor	31	MM227	Randolph Cyclodialysis Cannula, 19ga.
14	S0083	Graefe Muscle Hook #2	32	PV-S	PVA Sponges box of 50
15	S0079	Gass Retinal Detachment Hook	33	S0022	Bipolar Cord
16	S0037	Castroviejo Caliper, straight	34	S0023	Bipolar Forceps, angled blunt
17	S0080	Gass Scleral Marker	35	S0159	Squeeze Handle Forceps, serrated, straight
18	S0150	Schocket Scleral Depressor	36	S0160	Squeeze Handle Scissors, horizontal, straight
			37	S0161	Sterilizing Case, plastic (<i>not illustrated</i>)


CORNEAL TRANSPLANT SURGERY



AOLSS6- CORNEAL TRANSPLANT SURGERY

1	S0025	Bishop-Harmon A/C Irrigator	15	K20-2108	Barron Vacuum Punch, 8.0mm
2	S0103	Jones Towel Clamp	16	K21-2110	Barron Vacuum Punch, 8.5mm (not shown)
3	S0077	Flieringa Rings, (set of 8)	17	K5-1560	Polack Double Corneal Fixation Forceps
4	S0131	McNeill Global Fixation Ring	18	S0040	Castroviejo Forceps, 0.12mm
5	S0019	Barraquer Wire Speculum, large	19	S0041	Castroviejo Forceps, 0.3mm
6	S0119	Lieberman Speculum, solid blade	20	S0088	Harms-Colibri Forceps, 0.12mm
7	S0089	Hartmann Mosquito Forceps, straight (2)	21	S0169	Tennant Tying Forceps, straight
8	S0152	Serrefine, small, straight (2)	22	S0168	Tennant Tying Forceps, curved
9	S0175	Troutman-Katzin Scissors, small, right	23	S0143	Paton Spatula & Spoon
10	S0174	Troutman-Katzin Scissors, small, left	24	S0140	Needle Holder, curved w/lock
11	S0186	Westcott Tenotomy Scissors	25	S0043	Castroviejo Needle Holder, curved w/lock
12	S0098	Jaffe Stitch Scissors	26	S0003	Air Injection Cannula, 30ga.
13	K20-2056	Barron Vacuum Trepphine, 7.5mm	27	MM051-26	McIntyre A/C Cannula, 26ga., angled
14	K205258	Barron Vacuum Trepphine, 8.0mm (not shown)	28	S0161	Sterilizing Case, plastic (not illustrated)
			29	PV-S	PVA Sponges box of 50 (not illustrated)


STRABISMUS SURGERY



AOLSS7 - STRABISMUS SURGERY as used by Eugene M. Helveston, M.D

1 S0119	Lieberman Speculum, solid blade, small	15 S0037	Castroviejo Caliper, straight
2 S0216	Lieberman Speculum, solid blade, large	16 S0186	Westcott Tenotomy Scissors
3 S0089	Hartmann Mosquito Forceps, curved (2)	17 S0134	Westcott Stitch Scissors
4 S0152	Serrefine, small, straight (2)	18 S0029	Bonn Forceps, 0.12mm
5 S0212	Helveston Retractor, 7mm wide	19 S0217	Moody Fixation Forceps, left
6 S0213	Helveston Retractor, 9mm wide (not shown)	20 S0218	Moody Fixation Forceps, right
7 S0092	Helveston Retractor, 11mm, concave (not shown)	21 S0113	Lester Fixation Forceps, 2x3 (3)
8 S0215	Helveston Teaser Hook, 6mm (3)	22 S0093	Helveston Tying Forceps (2)
9 S0207	Helveston Finder Hook, large	23 S0043	Castroviejo Needle Holder, straight, w/lock
10 S0208	Helveston Finder Hook, small	24 S0206	Barraquer Needle Holder, curved, w/lock
11 S0211	Helveston Muscle Hook, 8mm	25 S0025	Bishop-Harmon A/C Irrigator
12 S0209	Helveston Muscle Hook, 10mm	26 S0051	Chuck Handle
13 S0210	Helveston Muscle Hook, 12mm	27 S0219	Performance Miniature Blade #69
14 S0214	Helveston Scleral Ruler	28 S0161	Sterilizing Case, plastic (not illustrated)

GLAUCOMA SURGERY



AOLSS8 - GLAUCOMA SURGERY

- | | | | |
|-----------|-------------------------------------|--------------|-------------------------------------|
| 1. S0025 | Bishop-Harmon A/C Irrigator | 16. S0060 | Diamond Step Knife, 10-facet |
| 2. S0103 | Jones Towel Clamp | 17. S0051 | Chuck Hand e (2) |
| 3. S0089 | Hartmann Mosquito Forceps, straight | 18. S0144 | Performance Lamellar Blade #66 (bx) |
| 4. S0151 | Schott Eye Speculum, adjustable | 19. MM-1500 | Performance Super Sharp Blade (bx) |
| 5. S0019 | Barraquer Wire Speculum, large | 20. I437 | Marking Pen (pk/10) |
| 6. S0070 | Eye Scissors, straight | 21. PVS-S | K-Sponge Spear, sterile (bx) |
| 7. S0186 | Westcott Tenotomy Scissors | 22. S0122 | Luntz-Dodick Trabeculectomy Punch |
| 8. S0016 | Barraquer Iris Scissors | 23. S0027 | Blumenthal Conjunctiva Dissector |
| 9. S0178 | Katena-Vannas Scissors, 7mm, curved | 24. S0018 | Barraquer Needle Holder, curved |
| 10. S0134 | McPherson-Westcott Stitch Scissors | 25. MM051-27 | Air Injection Cannula, 27ga. |
| 11. S0071 | Fechtner Conjunctiva Forceps (2) | 26. S0087 | Harms Trabeculectomy Probe, right |
| 12. S0053 | Colibri Forceps, 0.12mm teeth | 27. S0086 | Harms Trabeculectomy Probe, left |
| 13. S0026 | Bishop-Harmon Forceps | 28. S0161 | Sterilizing Case (not illustrated) |
| 14. S0072 | Fechtner Trabeculectomy Marker | | |
| 15. S0036 | Castroviejo Caliper | | |

ENUCLEATION SURGERY



AOLSS9- ENUCLEATION SURGERY

- | | | |
|----|--------------|--|
| 1. | S0115 | Lester-Burch Eye Speculum |
| 2. | S0184 | Wells Enucleation Spoon |
| 3. | S0083 | Graefe Muscle Hook #2 |
| 4. | S0165 | Steavens Scissors, straight |
| 5. | S0064 | Enucleation Scissors, strong curve |
| 6. | S0170 | Tissue Forceps |
| 7. | S0084 | Halsted Mosquito Forceps, curved |
| 8. | S0161 | Sterilizing Case, plastic (<i>not illustrated</i>) |


LASIK SURGERY



AOLSS10 - LASIK SURGERY

- | | |
|-----------|--|
| 1. S0118 | Lieberman Speculum, K-wire |
| 2. S0124 | Buratto Flap Protector |
| 3. S0124 | Mannis-Buratto Flap Protector |
| 4. S0094 | Hersh LASIK Spatula |
| 5. S0105 | Katzen LASIK Flap "Unzipper" |
| 6. S0082 | Gimbe Stabilization Ring |
| 7. S0050 | Chayet Corneal Marker |
| 8. S0033 | Buratto Flap Forceps |
| 9. S0158 | Spatulated Irrigation Cannu a, 26ga. |
| 10. S0034 | Buratto Irrigation Cannula, 3-ports |
| 11. S0102 | John IntraLasik Cannu a, sterie (pk/10) |
| 12. PV-S | PVA-Sponge pack of 50 |
| 13. I437 | Marking Pen, standard (pk/10) |
| 14. S0117 | Lieberman Aspiration Speculum |
| 15. S0161 | Sterilizing Case, plastic (<i>not illustrated</i>) |
- * optional

LASEK SURGERY



AOLSSI I - LASEK SURGERY

- | | | |
|----|--------------|--|
| 1. | S0117 | Lieberman Aspirating Speculum |
| 2. | PV-S | PVA-Sponge pack of 50 |
| 3. | S0004 | Alcohol Cannula |
| 4. | S0172 | TrephineS Well 8/8.5mm |
| 5. | S0173 | Trephine & Wei 9/9.5mm |
| 6. | S0139 | Micro-Hoe & Epi-Peeler |
| 7. | S0076 | Flap Repositor |
| 8. | S0161 | Sterilizing Case, plastic (<i>not illustrated</i>) |

LID SURGERY



AOLSS I2- LID SURGERY

- | | | | |
|------------|-------------------------------------|-------------|---|
| 1. S0025 | Bishop-Harmon A/C Irrigator | 17. S0056 | Desmarres Lid Retractor #2 |
| 2. S0085 | Halsted Mosquito Forceps, straight | 18. S0057 | Desmarres Lid Retractor #3 |
| 3. S0084 | Halsted Mosquito Forceps, curved (2 | 19. S0083 | Graefe Muscle Hook #2 |
| 4. S0103 | Jones Towel Clamp, 2F | 20. S0075 | Fixation Hook, double |
| 5. S0165 | Stevens Scissors, straight | 21. S0028 | Bonaccolto Conjunctiva Forceps |
| 6. S0070 | Eye Scissors, straight | 22. S0026 | Bishop-Harmon Forceps |
| 7. S0051 | Chuck Han die | 23. S0041 | Castroviejo Forceps, 0.3mm |
| 8. S0147 | Performance Miniature B ade #67 (bx | 24. S0093 | Helveston Tying Forceps |
| 9. S0186 | Westcott Tenotomy Scissors | 25. S0062 | Dowries Lid Clamp |
| 10. S0134 | Westcott Stitch Scissors | 26. KI-8481 | Scleral Shield, small, pair |
| 11. S0021 | Berke Ptosis Forceps | 27. KI-8482 | Scleral Shield, large, pair |
| 12. S0018 | Barraquer Needle Holder, curved | 28. S0023 | Bipolar Forceps, angled, blunt |
| 13. K-8520 | Jaeger Lid Plate, stain ess steel | 29. S0022 | Bipolar Cord |
| 14. S0042 | Castroviejo Needle Holder, curved | 30. S0090 | Helveston Retractor, concave, 11 mm |
| 15. S0037 | Castroviejo Caliper, straight | 31. S0109 | Knapp Retractor, 4-prong |
| 16. S0055 | Desmarres Lid Retractor #1 | 32. I437 | Marking Pen, standard (pk/10) |
| | | 33. S0161 | Sterilizing Case, plastic (not illustrated) |

CHALAZION SURGERY



AOLSS13- CHALAZION SURGERY

1. **KI-8520** Jaeger Lid Plate, stainless steel
2. **S0135** Meyerhoefer Curette #1
3. **S0136** Meyerhoefer Curette #2
4. **S0137** Meyerhoefer Curette #3
5. **S0138** Meyerhoefer Curette #4
6. **S0067** Eye Scissors, curved
7. **S0123** Manhattan Eye & Ear Suture Forceps
8. **S0095** Hunt Chalazion Forceps
9. **S0078** Francis Chalazion Forceps
10. **S0042** Castroviejo Needle Holder, curved
11. **S0147** Performance Miniature Blade #67 (bx)
12. **S0051** Chuck Handle
13. **S0161** Sterilizing Case, plastic (*not illustrated*)

PTERYGIUM SET



AOLSSI 4- PTERYGIUM SET

- | | | |
|-----|----------------|-----------------------------------|
| 1. | S0145 | Miniature B ade #64, disposable |
| 2. | S0051 | Chuck Handle |
| 3. | S0112 | Lancaster Eye Speculum |
| 4. | S0133 | McPherson Tying Forceps, straight |
| 5. | S0040 | Castroviejo Forceps, 0.12mm |
| 6. | S0081 | Gi I Corneal Knife |
| 7. | S0042 | Castroviejo Needle Holder, curved |
| 8. | S0142 | Katena Ophthalmic Burr |
| 9. | K2-4925 | Diamond Burr, 5.0mm diameter ba |
| 10. | S0070 | Eye Scissors, straight, standard |
| 11. | S0068 | Eye Scissors, straight, smal |
| 12. | S0028 | Bonaccolto Conjunctiva Forceps |
| 13. | S0113 | Lester Fixation Forcep |

VETERINARY OPHTHALMIC SURGERY SET



AOLSS15-VETERINARY OPHTHALMIC SURGERY SET

- | | | | |
|------------------|--|------------------|---|
| 1. S0020 | Barraquer Wire Speculum, small | 15. S0097 | Intra-Ocular Scissors, straight |
| 2. S0019 | Barraquer Wire Speculum, arge | 16. S0181 | Intra-Ocular Scissors, angled 30* |
| 3. S0030 | Iris Forceps | 17. S0001 | Adson Tissue Forceps |
| 4. S0052 | Colibri Forceps, 1x2 teeth, 0.4mm | 18. S0044 | Castroviejo Needle Holder, delicate, curved |
| 5. S0088 | Harms-Colibri Forceps, 0.12mm | 19. S0045 | Castroviejo Needle Holder, heavy, curved |
| 6. S0039 | Castroviejo Corneal Scissors, small, right | 20. S0157 | Sinsky Lens Hook, angled |
| 7. S0038 | Castroviejo Corneal Scissors, small, left | 21. S0099 | Jaffe-Bechert Rotator, 9mm |
| 8. S0186 | Westcott Tenotomy Scissor | 22. S0051 | Chuck Handle (2) |
| 9. S0177 | Utrata Capsulorrhexis Forcep | 23. S0015 | Bard-Parker Handle |
| 10. S0026 | Bishop-Harmon Forceps | 24. S0103 | Jones Towel C amp (4) |
| 11. S0063 | Dressing Forceps, delicate, straight | 25. S0085 | Halsted Mosquito Forceps, straight (2) |
| 12. S0171 | Tissue Forceps, delicate, straight | 26. S0084 | Halsted Mosquito Forceps, curved (2) |
| 13. S0049 | Chang Phaco Chopper | 27. S0141 | Utility Scissors |
| 14. S0182 | Warren Canine Phaco Chopper | 28. Mm080 | Chang Hydrodissection Cannula |

Head Office
Akriti Ophthalmic LLC,
330 N BRAND BLVD STE 700
GLENDALE, CA 91203
USA
+1 (864) 624 6898
sales@akritiophthalmic.com
www.akritiophthalmic.com

Asia
Akriti Ophthalmic Pvt Ltd,
M34, M35, M36 Medical Devices Park,
Sultanpur Pathancheru,
Hyderabad - 502319,
Telangana State, India.
+91 995 954 1995
sales@akriti.co.in
www.akriti.co.in

Middle East & Africa
Akriti Medical Equipments Trading Co L.L.C
The Prism Tower,
Flat No 3201, Business bay, Dubai,
United Arab Emirates
+971 547 917 448
sales@akritigroup.com
www.akritigroup.com